



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05542 - Goat Anti-ARL1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** ARL1, ADP-ribosylation factor-like 1, ARFL1

**Official Symbol:** ARL1

**Accession Number(s):** NP\_001168

**Human GeneID(s):** [400](#)

**Non-Human GeneID(s):** 66941 (mouse)

### Immunogen

Peptide with sequence C-EAMEWLVELKSRQ, from the C Terminus of the protein sequence according to NP\_001168.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Western blot:** Western blot: Approx 20kDa band seen in Human Pancreas lysate [Predicted MW of 20kDa according to NP\_001168]. Recommended for use at 1-3µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat