

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB07788 - Goat Anti-ASF1A / HSPC146 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ASF1A, HSPC146, ASF1 anti-silencing function 1 homolog A (S. cerevisiae), RP3-329L24.1, CGI-98, CIA, DKFZP547E2110, ASF1 anti-silencing function

1 homolog A, CCG1-interacting factor A, anti-silencing function 1A

Official Symbol: ASF1A

Accession Number(s): NP_054753.1

Human GeneID(s): 25842

Non-Human GenelD(s): 66403 (mouse), 294408 (rat)

Immunogen

Peptide with sequence C-EDNTEKLEDAE, from the internal region of the protein sequence according to NP_054753.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Approx. 23kDa band observed in lysates of cell line Jurkat (calculated MW of 23.0kDa according to NP_054753.1). Recommended concentration: 0.01-0.03µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog



EB07788 ($0.01\mu g/ml$) staining of Jurkat cell lysate ($35\mu g$ protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.