

International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07635 - Goat Anti-FRAG1 (mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Atad5, ATPase family, AAA domain containing 5, C17orf41, chromosome 17 open reading frame 41, FLJ12735, FRAG1, chromosome fragility associated gene 1, hypothetical protein LOC237877, RP23-328O2.1, C130052G03Rik

Official Symbol: Atad5

Accession Number(s): NP_001025027.1

Non-Human GeneID(s): 237877 (mouse), 303348 (rat)

Immunogen

Peptide with sequence C-QSKDQPLRKSQK, from the internal region of the protein sequence according to NP_001025027.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary experiments gave an approx 60kDa band in lysates of cell lines NIH3T3 and NSO after 0.01µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 204kDa according to NP_001025027.1. The 60kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse, Rat