



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08862 - Goat Anti-ATG5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ATG5, ATG5 autophagy related 5 homolog (*S. cerevisiae*), APG5, APG5L, ASP, hAPG5, APG5 (autophagy 5, *S. cerevisiae*)-like, APG5 autophagy 5-like, ATG5 autophagy related 5-like, OTTHUMP00000040507, apoptosis specific protein

Official Symbol: ATG5

Accession Number(s): NP_004840.1; NP_001273036.1; NP_001273037.1

Human GeneID(s): [9474](#)

Non-Human GeneID(s): 11793 (mouse)

Immunogen

Peptide with sequence C-SQVINEMQKKDHKQ, from the internal region of the protein sequence according to NP_004840.1; NP_001273036.1; NP_001273037.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary testing showed a band at approx. 28-30kDa in lysate of cell line MDA-MB-231, after 1µg/ml antibody staining (calculated Mwt. Of 32.4 kDa according to NP_004840.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Cow