



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05440 - Goat Anti-ATP11B Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ATP11B, ATPase, Class VI, type 11B, ATP1F, ATP1R, KIAA0956, DKFZP434J238, DKFZP434N1615

Official Symbol: ATP11B

Human GeneID(s): [23200](#)

Immunogen

Peptide with sequence DRSILTLSTMDSSSTC, from the C Terminus of the protein sequence according to .

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 85-90kDa band seen in 293 lysate. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human