

### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05440 - Goat Anti-ATP11B Antibody

Size: 100µg specific antibody in 200µl



### **Target Protein**

Principal Names: ATP11B, ATPase, Class VI, type 11B, ATPIF, ATPIR, KIAA0956,

DKFZP434J238, DKFZP434N1615

Official Symbol: ATP11B Human GeneID(s): 23200

### Immunogen

Peptide with sequence DRSILTLSTMDSSTC, from the C Terminus of the protein sequence according to .

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

# **Applications Tested**

**Western blot:** Western blot: Approx 85-90kDa band seen in 293 lysate. Recommended for use at 1-3μg/ml.

# **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human