

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05402 - Goat Anti-BAP29 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BAP29, B-cell receptor-associated protein BAP29

Official Symbol: BCAP29

Accession Number(s): NP_061332

Human GenelD(s): 55973

Immunogen

Peptide with sequence CELQDRLERGNKKRL, from the C Terminus of the protein sequence according to NP_061332.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 25kDa band seen in Mouse Spleen lysate. Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human