



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB12083 - Goat Anti-BATF Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BATF, basic leucine zipper transcription factor, ATF-like, B-ATF, BATF1, SFA-2, SFA2, B-cell-activating transcription factor, SF-HT-activated gene 2 protein, activating transcription factor B, basic leucine zipper transcriptional factor ATF-like

Official Symbol: BATF

Accession Number(s): NP_006390.1

Human GeneID(s): [10538](#)

Non-Human GeneID(s): 53314 (mouse), 299206 (rat)

Immunogen

Peptide with sequence C-EVVYSAHAFHQPHVS, from the C Terminus of the protein sequence according to NP_006390.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Preliminary experiments gave an approx 25kDa band in Human Spleen lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 14.1kDa according to NP_006390.1. The 25kDa band was successfully blocked by incubation with the immunizing peptide.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow