



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05789 - Goat Anti-BAZ2A / TIP5 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BAZ2A, bromodomain adjacent to zinc finger domain, 2A, TIP5, KIAA0314, TTF-I interacting peptide 5, DKFZp781B109, FLJ13768, FLJ13780, FLJ45876, WALp3

Official Symbol: BAZ2A

Accession Number(s): NP_038477.2

Human GeneID(s): [11176](#)

Non-Human GeneID(s): 116848 (mouse)

Immunogen

Peptide with sequence C-RWEEFYQGKQANL, from the C Terminus of the protein sequence according to NP_038477.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:2000.

Western blot: Western Blot: No signal obtained yet but low background observed in 293, U937, NIH-3T3, Human Brain and Human Heart lysates at upto 2µg/ml.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow