

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB10601 - Goat Anti-Bbs1 (mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Bbs1, Bardet-Biedl syndrome 1, Al451249, D19Ertd609e, Bardet-Biedl

syndrome 1, Bardet-Biedl syndrome 1 homolog

Official Symbol: Bbs1

Accession Number(s): NP_001028300.1

Human GeneID(s): 582

Non-Human GenelD(s): 52028 (mouse), 309156 (rat)

Immunogen

Peptide with sequence C-NGSIYILRRDSKHPK, from the internal region of the protein sequence according to NP_001028300.1.

Please note the <u>peptide</u> is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog