



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08340 - Goat Anti-BDH2 / DHRS6 (aa 92 to 103) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BDH2, 3-hydroxybutyrate dehydrogenase, type 2, DHRS6, EFA6R, FLJ13261, PRO20933, UCPA-OR, UNQ6308, dehydrogenase/reductase (SDR family) member 6, oxidoreductase UCPA

Official Symbol: BDH2

Accession Number(s): NP_064524.3

Human GeneID(s): [56898](#)

Non-Human GeneID(s): 69772 (mouse), 295458 (rat)

Immunogen

Peptide with sequence DCEEKDWDFSMN, from the internal region of the protein sequence according to NP_064524.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Preliminary experiments in Human Brain (Cerebellum), Liver and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Cow