



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05760 - Goat Anti-BIRC6 / Apollon Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BIRC6, baculoviral IAP repeat-containing 6 (apollon), KIAA1289, baculoviral IAP repeat-containing 6, BRUCE, APOLLON, FLJ13726, FLJ13786, KIAA1289, ubiquitin-conjugating BIR-domain enzyme apollon

Official Symbol: BIRC6

Accession Number(s): NP_057336.3

Human GeneID(s): [57448](#)

Immunogen

Peptide with sequence C-PSSSKELPSDFQL, from the C Terminus of the protein sequence according to NP_057336.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.


Western blot: Preliminary experiments gave no signal but low background in Human Liver and Kidney lysates at up to 1µg/ml. However this used our routine western blotting protocol which we would not expect to reliably detect proteins as large as the predicted size of approx. 555KDa according to NP_057336.

IHC: Paraffin embedded Human Lung. Recommended concentration: 3.75µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB05760 (3.75µg/ml) staining of paraffin embedded Human Lung. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.