



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB09000 - Goat Anti-CCL3L1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CCL3L1, chemokine (C-C motif) ligand 3-like 1, 464.2, D17S1718, G0S19-2, LD78, LD78BETA, MGC104178, MGC12815, SCYA3L, SCYA3L1, small inducible cytokine A3-like 1

**Official Symbol:** CCL3L1

**Accession Number(s):** NP\_066286.1

**Human GeneID(s):** [6349](#)

### Immunogen

Peptide with sequence CNQVLSAPLAADT , from the internal region of the protein sequence according to NP\_066286.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:8000.

**Western blot:** CCL3L1 expression requires activation.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human