

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB08119 - Goat Anti-CCRN4L / Nocturnin Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** CCRN4L, CCR4 carbon catabolite repression 4-like (S. cerevisiae), CCR4L, MGC142054, MGC142060, MGC4120817, MGC78549, NOC, CCR4 carbon catabolite repression 4-like, CCR4-like (carbon catabolite repression 4, S.cerevisiae), nocturnin

Official Symbol: CCRN4L

Accession Number(s): NP\_036250.2

Human GenelD(s): 25819

### **Immunogen**

Peptide with sequence HLSLVCDFSFTEESD, from the C Terminus of the protein sequence according to NP\_036250.2.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:4000.

**Western blot:** Western Blot: Approx 37kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 48.2kDa according to NP\_036250.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Garbarino-Pico et al, jRNA. 2007 May;13(5):745-55; PMID: 17400819). Recommended concentration: 0.3-1µg/ml.

## **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human