



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10965 - Goat Anti-CD127 / IL7R (mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: IL7R, interleukin 7 receptor, CD127, IL-7R-alpha, IL7RA, ILRA, CD127 antigen, IL-7 receptor subunit alpha, IL-7R subunit alpha, interleukin 7 receptor alpha chain, interleukin 7 receptor isoform H5-6, interleukin-7 receptor subunit alpha, MGC10755

Official Symbol: IL7R

Accession Number(s): NP_032398.3

Non-Human GeneID(s): 16197 (mouse), 294797 (rat)

Immunogen

Peptide with sequence DDHSFWCHSQLE, from the internal region (near N Terminus) of the protein sequence according to NP_032398.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary testing showed a band at approx 50kDa in NIH3T3 cell lysate, and 55kDa in Mouse Duodenum lysate after 1-2µg/ml antibody staining (calculated MW of 51.6kDa according to Mouse NP_032398.3). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Mouse, Rat