



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08288 - Goat Anti-CDKAL1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CDK5 regulatory subunit associated protein 1-like 1, FLJ20342, FLJ46705, MGC75469

Official Symbol: CDKAL1

Accession Number(s): NP_060244.2

Human GeneID(s): [54901](#)

Non-Human GeneID(s): 68916 (mouse), 361243 (rat)

Immunogen

Peptide with sequence C-QVKKQRTKDLSR, from the internal region of the protein sequence according to NP_060244.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: This antibody has been described in a publication and it has been successfully used on Human islets and pancreatic insulinoma INS-1 and MIN6 cells. Preliminary experiments gave an approx. 39kDa band in Human Brain (Cerebellum) lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 65.1kDa according to NP_060244.2. The 39kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog

Specific Reference

Brambillasca S, Altkrueger A, Colombo SF, Friederich A, Eickelmann P, Mark M, Borgese N, Solimena M.

CDK5 regulatory subunit-associated protein 1-like 1 (CDKAL1) is a tail-anchored protein in the endoplasmic reticulum (ER) of insulinoma cells.

J Biol Chem. 2012 Dec 7;287(50):41808-19.

PMID: 23048041