

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05621 - Goat Anti-Centromere protein E Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CENPE, CENP-E, centromere protein E (312kD), Centromere autoantigen E (312kD), centromere protein E, 312kDa, KIF10, centromere protein E,

kinesin family member 10 **Official Symbol:** CENPE

Accession Number(s): NP_001804.2

Human GeneID(s): 1062

Immunogen

Peptide with sequence HASSGKDVPECKTQ, from the C Terminus of the protein sequence according to NP_001804.2.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the predicted size of approx. 316kDa according to NP_001804. Therefore we cannot recommend an optimal concentration and the product is Investigative Grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Rat, Dog, Cow