

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08477 - Goat Anti-CHRNA4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CHRNA4, cholinergic receptor, nicotinic, alpha 4, BFNC, EBN, EBN1, NACRA4, cholinergic receptor, nicotinic, alpha 4 subunit, cholinergic receptor, nicotinic, alpha polypeptide 4, neuronal acetylcholine receptor alpha 4 subunit, neuronal nicotinic

acetylcholine receptor alpha-4 subunit

Official Symbol: CHRNA4

Accession Number(s): NP_000735.1

Human GeneID(s): 1137

Immunogen

Peptide with sequence C-RASSHVETRAHAE, from the internal region (near the N Terminus) of the protein sequence according to NP_000735.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Preliminary experiments in Human Brain (Cerebral Cortex, Cerebellum and Frontal Cortex) lysates gave no specific signal but low background (at antibody concentration up to 1μg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human