



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB11018 - Goat Anti-Col11a2 (aa256-268) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Col11a2, collagen, type XI, alpha 2, RP23-222K15.2, collagen alpha-2(XI) chain, procollagen, type XI, alpha 2, collagen alpha-2(XI) chain, collagen, type XI, alpha 2, procollagen, type XI, alpha 2, RP23-222K15.2, Col11a2

Official Symbol: Col11a2

Accession Number(s): NP_034056.1

Non-Human GeneID(s): 12815 (mouse)

Immunogen

Peptide with sequence C-HKPQSQEPQKQYP, from the internal region of the protein sequence according to NP_034056.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Species Reactivity

Tested:

Expected from sequence similarity: Mouse