

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09671 - Goat Anti-Connexin 30 / GJB6 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GJB6, gap junction protein, beta 6, 30kDa, CX30, DFNA3, ED2, EDH, HED, OTTHUMP0000018096, OTTHUMP00000176870, OTTHUMP00000176871, OTTHUMP00000176872, connexin 30, ectodermal dysplasia 2, hidrotic (Clouston syndrome), gap junction protein, beta 6, gap junction protein, beta 6 (connexin 30)

Official Symbol: GJB6

Accession Number(s): NP_006774.2

Human GeneID(s): [10804](#)

Non-Human GeneID(s): 14623 (mouse), 84403 (rat)

Important Comments: Reported variants represent identical protein: NP_006774.2, NP_001103691.1, NP_001103689.1, NP_001103690.1

Immunogen

Peptide with sequence C-RSKRAQTQKNHPNH, from the internal region of the protein sequence according to NP_006774.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave an approx 14kDa band in Human Brain (Cerebellum) and Rat Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 30.4kDa according to NP_006774.2. The 14kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Cow