

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09672 - Goat Anti-Connexin 50 / GJA8 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GJA8, gap junction protein, alpha 8, 50kDa, CAE, CAE1, CX50, CZP1, MP70, OTTHUMP00000025757, connexin 50, gap junction membrane channel protein alpha-8, connexin-50, gap junction alpha-8 protein, lens fiber protein MP70, lens intrinsic membrane protein MP70

Official Symbol: GJA8

Accession Number(s): NP_005258.2

Human GeneID(s): [2703](#)

Non-Human GeneID(s): 14616 (mouse), 29601 (rat)

Immunogen

Peptide with sequence TPSLCPELTDDAR, from the internal region of the protein sequence according to NP_005258.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Preliminary experiments in Human Brain (Spinal Cord), Mouse and Rat Eye lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow