



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05116 - Goat Anti-CRK isoform II Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: CRK, v-crk sarcoma virus CT10 oncogene homolog (avian), CRKII, v-crk avian sarcoma virus CT10 oncogene homolog, avian sarcoma virus CT10 (v-crk) oncogene homolog

Official Symbol: CRK

Accession Number(s): NM_016823

Human GeneID(s): [1398](#)

Important Comments: This antibody is expected to recognise the second human isoform only.

Immunogen

Peptide with sequence C-RLLDQQNPDEDFS, from the C Terminus of the protein sequence according to NM_016823.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 35kDa band seen in Human Brain lysate.

Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse