



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05118 - Goat Anti-CSK Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CSK, c-src tyrosine kinase, MGC117393

**Official Symbol:** CSK

**Accession Number(s):** NP\_004374.1

**Human GeneID(s):** [1445](#)

**Non-Human GeneID(s):** 12988 (mouse), 315707 (rat)

**Important Comments:** Reported variants represent identical protein: NP\_004374.1; NP\_001120662.1.

### Immunogen

Peptide with sequence C-EQLEHIKTHELHL, from the C Terminus of the protein sequence according to NP\_004374.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested


**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Approx 48kDa band observed in lysates of Jurkat, K562 and MOLT4 (calculated MW of 50.7kDa according to NP\_004374.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Cow



EB05118 (0.1µg/ml) staining of Jurkat lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.