



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06401 - Goat Anti-CYLD (C terminus) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** CYLD, EAC, CDMT, CYLD1, CYLDI, FLJ31664, FLJ78684, HSPC057, KIAA0849, USPL2, deubiquitinating enzyme CYLD, ubiquitin carboxyl-terminal hydrolase CYLD, ubiquitin specific peptidase like 2, ubiquitin thiolesterase CYLD, ubiquitin-specific-processing protease CYLD, cylindromatosis (turban tumor syndrome), MFT, MFT1, SBS, TEM

**Official Symbol:** CYLD

**Accession Number(s):** NP\_056062.1; NP\_001035814.1; NP\_001035877.1

**Human GeneID(s):** [1540](#)

**Important Comments:** This antibody is expected to recognize both reported isoforms (NP\_056062.1 and NP\_001035814.1; NP\_001035877.1).

### Immunogen

Peptide with sequence CMYQSPTMSLYK, from the C Terminus of the protein sequence according to NP\_056062.1; NP\_001035814.1; NP\_001035877.1.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:16000.

**Western blot:** Preliminary experiments gave no signal but low background in human kidney and Hela lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Cow

### Background Reference

Regamey A, Hohl D, Liu JW, Roger T, Kogerman P, Toftgard R, Huber M.

The tumor suppressor CYLD interacts with TRIP and regulates negatively nuclear factor kappaB activation by tumor necrosis factor.

J Exp Med. 2003 Dec 15;198(12):1959-64.

PMID: 14676304