



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07489 - Goat Anti-DAP10 / HCST Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HCST, DAP10, hematopoietic cell signal transducer, KAP10, PIK3AP, DNAX-activation protein 10, phosphoinositide-3-kinase adaptor protein

Official Symbol: HCST

Accession Number(s): NP_055081.1; NP_001007470.1

Human GeneID(s): [10870](#)

Non-Human GeneID(s): 23900 (mouse), 474146 (rat)

Important Comments: This antibody is expected to recognise isoform 1 (NP_055081.1) but it may recognize isoform 2 lacking residue E81 (NP_001007470.1).

Immunogen

Peptide with sequence C-RSPAQEDGKVVYIN, from the C Terminus of the protein sequence according to NP_055081.1; NP_001007470.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat