

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB09922 - Goat Anti-DDA1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: C19orf58, cross-immune reaction antigen PCIA1, DET1 and DDB1

associated 1, MGC2594, PCIA1, DDA1

Official Symbol: DDA1

Accession Number(s): NP_076955.1

Human GenelD(s): 79016

Non-Human GenelD(s): 66498 (mouse)

Immunogen

Peptide with sequence C-PTREYPSEQIIVTE, from the internal region of the protein sequence according to NP_076955.1.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave an approx 25kDa band in Mouse Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 11.8kDa according to NP_076955.1. The 25kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Dog

Background Reference

Olma MH, Roy M, Le Bihan T, Sumara I, Maerki S, Larsen B, Quadroni M, Peter M, Tyers M, Pintard L,

An interaction network of the mammalian COP9 signalosome identifies Dda1 as a core subunit of multiple Cul4-based E3 ligases.

Journal of cell science 2009 Apr 122 (Pt 7): 1035-44.

PMID: 19295130