

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

#### www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB09461 - Goat Anti-Decorin Antibody

Size: 100µg specific antibody in 200µl



### **Target Protein**

**Principal Names:** DCN, decorin, CSCD, DSPG2, PG40, PGII, PGS2, SLRR1B, bone proteoglycan II, decorin proteoglycan, dermatan sulphate proteoglycans II, proteoglycan core protein, small leucine-rich protein 1B

Official Symbol: DCN

Accession Number(s): NP\_001911.1; NP\_598011.1; NP\_598012.1;

Human GeneID(s): 1634

**Important Comments:** This antibody is expected to recognize the three reported isoforms a, b and c (NP\_001911.1; NP\_598011.1; NP\_598012.1 resp.). Reported variants NP\_598010.1 and NP\_001911.1 represent identical protein.

# **Immunogen**

Peptide with sequence C-KISRVDAASLKGLNN, from the internal region of the protein sequence according to NP\_001911.1; NP\_598011.1; NP\_598012.1;.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:8000.

**Western blot:** Approx 38kDa band observed in Human Kidney lysates (calculated MW of 39.7kDa according to NP\_001911.1). Recommended concentration: 0.5-2μg/ml. Primary incubation was 1 hour.

IHC: Paraffin embedded Human Breast. Recommended concentration: 3.75µg/ml.

### **Species Reactivity**


Tested: Human

Expected from sequence similarity: Human, Cow, Pig, Cavia



uining of Human Kidnov lysata (25ug protain in Pl

EB09461 (0.5 $\mu$ g/ml) staining of Human Kidney lysate (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.



EB09461 (3.75 $\mu$ g/ml) staining of paraffin embedded Human Breast. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.