

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB11108 - Goat Anti-DEFB104A / DEFB104B (aa37-50) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: BD-4, beta-defensin 104, beta-defensin 4, DEFB104, DEFB104B, DEFB4, DEFB-4, defensin, beta 104A, defensin, beta 4, hBD-4, MGC118942,

MGC118944, MGC118945, DEFB104A

Official Symbol: DEFB104A

Accession Number(s): NP_525128.2; NP_001035792.1

Human GenelD(s): 140596, 503618

Immunogen

Peptide with sequence CRKKCRSQEYRIGR, from the internal region of the protein sequence according to NP_525128.2; NP_001035792.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins smaller than the calculated size of 8kDa according to the precursor NP_525128.2;. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product.

Species Reactivity

Tested:

Expected from sequence similarity: Human