



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05211 - Goat Anti-DOK4 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** DOK4, downstream of tyrosine kinase 4, FLJ10488

**Official Symbol:** DOK4

**Accession Number(s):** NP\_060580

**Human GeneID(s):** [55715](#)

### Immunogen

Peptide with sequence C-SQGDSSEAKTPSQ, from the C Terminus of the protein sequence according to NP\_060580.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Western blot:** Western blot: Approx 45-50kDa band seen in Mouse Kidney extracts. Recommended for use at 0.5-2µg/ml. Overnight incubation is recommended with this antibody.

### Species Reactivity

**Tested:** Mouse

**Expected from sequence similarity:** Human