



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05186 - Goat Anti-DRAK2 / STK17B Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** STK17B, serine/threonine kinase 17b (apoptosis-inducing), DRAK2, death-associated protein kinase-related 2

**Official Symbol:** STK17B

**Accession Number(s):** NP\_004217.1

**Human GeneID(s):** [9262](#)

### Immunogen

Peptide with sequence SLPNPHELVSDLLC, from the C Terminus of the protein sequence according to NP\_004217.1.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:8000.

**Western blot:** Western blot: Approx 45kDa band seen in Jurkat cell lysates.

Recommended for use at 2-4µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Pig, Cow