

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06156 - Goat Anti-DSS1 / SHFM1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SHFM1, split hand/foot malformation (ectrodactyly) type 1, DSS1, Deleted in split-hand/split-foot 1 region, ECD, SEM1, SHFD1, SHSF1, Shfdg1, candidate for split hand/foot malformation type 1, deleted in split-hand/foot 1, deleted in split-hand/split-foot 1, split hand/foot malformation type 1

Official Symbol: SEM1

Accession Number(s): NP_006295.1

Human GeneID(s): [7979](#)

Non-Human GeneID(s): 20422 (mouse)

Immunogen

Peptide with sequence C-RAELEKHGYKMETS, from the C Terminus of the protein sequence according to NP_006295.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: No signal obtained yet but a preliminary experiment showed low background in Human Liver extract at up to 1µg/ml. This used our standard western blotting protocol which we would not expect to detect proteins as small as the predicted size of approx. 8KDa according to NP_006295. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Immunocytochemistry: Anonymous customer found nuclear staining in HT1080.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Dog, Cow, Pig