



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06156 - Goat Anti-DSS1 / SHFM1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** SHFM1, split hand/foot malformation (ectrodactyly) type 1, DSS1, Deleted in split-hand/split-foot 1 region, ECD, SEM1, SHFD1, SHSF1, Shfdg1, candidate for split hand/foot malformation type 1, deleted in split-hand/foot 1, deleted in split-hand/split-foot 1, split hand/foot malformation type 1

**Official Symbol:** SEM1

**Accession Number(s):** NP\_006295.1

**Human GeneID(s):** [7979](#)

**Non-Human GeneID(s):** 20422 (mouse)

### Immunogen

Peptide with sequence C-RAELEKHGYKMETS, from the C Terminus of the protein sequence according to NP\_006295.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** No signal obtained yet but a preliminary experiment showed low background in Human Liver extract at up to 1µg/ml. This used our standard western blotting protocol which we would not expect to detect proteins as small as the predicted size of approx. 8KDa according to NP\_006295. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

**Immunocytochemistry:** Anonymous customer found nuclear staining in HT1080.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Dog, Cow, Pig