



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08652 - Goat Anti-DUBA2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: OTUD6A, OTU domain containing 6A, DUBA2, FLJ25831, HSHIN6, HIN-6 protease

Official Symbol: OTUD6A

Accession Number(s): NP_997203.1

Human GeneID(s): [139562](#)

Immunogen

Peptide with sequence RCRTASYMKKHVDE, from the internal region of the protein sequence according to NP_997203.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary experiments gave an approx 22-23kDa band in Human Pancreas, Human Ovary and Human Kidney lysates after 1µg/ml antibody staining.

Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 33.3kDa according to NP_997203.1. The 23kDa band was successfully blocked by incubation with the immunizing peptide.

Species Reactivity

Tested:

Expected from sequence similarity: Human