

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08378 - Goat Anti-Dudulin 2 / STEAP3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: STEAP3, STEAP family member 3, STMP3, TSAP6, dudlin-2, 1010001D01Rik, dudlin 2, dudulin 2, six transmembrane prostate protein 3, tumor suppressor pHyde

Official Symbol: STEAP3

Accession Number(s): NP_878919.2; NP_001008410.1; NP_060704.2

Human GeneID(s): [55240](#)

Important Comments: This antibody is expected to recognise all reported isoforms (NP_878919.2; NP_001008410.1 and NP_060704.2).

Immunogen

Peptide with sequence C-QPYVQESQNKFFK, from the internal region of the protein sequence according to NP_878919.2; NP_001008410.1; NP_060704.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Western Blot: Approx 50-55kDa band observed in Human Liver and Pancreas lysates (calculated MW of 55.7kDa according to NP_878919.2 and of 54.6kDa according to NP_001008410.1 and NP_060704.2). Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human