

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB11892 - Goat Anti-DUSP6 / MKP3 Antibody

Size: 100µg specific antibody in 200µl



# **Target Protein**

**Principal Names:** 

Official Symbol: DUSP6

Accession Number(s): NP\_001937.2; NP\_073143.2

Human GeneID(s): 1848

Non-Human GenelD(s): 67603 (mouse), 116663 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms

(NP\_001937.2; NP\_073143.2).

## **Immunogen**

Peptide with sequence C-PSNQNVYQVDSLQST, from the C Terminus of the protein sequence according to NP\_001937.2; NP\_073143.2.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**


Peptide ELISA: antibody detection limit dilution 1:128000.

**Western blot:** Approx 26kDa band observed in Mouse and Rat Heart lysates and 40kDa plus 26kDa bands observed in Human and Pig Heart lysates (calculated MW of 42.3kDa according to NP\_001937.2 and 26.5kDa according to NP\_073143.2). Recommended concentration: 0.3-1µg/ml.

### **Species Reactivity**


Tested: Human, Mouse, Rat, Pig

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB11892 (0.3μg/ml) staining of Human (A), (Mouse) and Rat (C) Heart lysates (35μg protein in RIPA buffer).

Primary incubation was 1 hour. Detected by chemiluminescence.



EB11892 (0.3μg/ml) staining of Pig Heart lysate (35μg protein in RIPA buffer). Primary incubation was 1 hour.

Detected by chemiluminescence.