

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07416 - Goat Anti-EDG2 / LPAR1 Antibody

Size: 100µg specific antibody in 200µl

Target Protein

Principal Names: EDG2, endothelial differentiation, lysophosphatidic acid

G-protein-coupled receptor, 2, Gpcr26, LPA1, LPAR1, Mrec1.3, edg-2, rec.1.3, vzg-1, ventricular zone gene 1

Official Symbol: EDG2

Accession Number(s): NP_001392.2; NP_476500.1

Human GeneID(s): [1902](#)

Non-Human GeneID(s): 14745 (mouse), 116744 (rat)

Important Comments: Both reported variants represent identical protein (NP_001392.2 and NP_476500.1).

Immunogen

Peptide with sequence C-SRHSSGPRRNRDT, from the internal region of the protein sequence according to NP_001392.2; NP_476500.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Preliminary experiments in lysates of Human Brain, Mouse Brain and Rat Brain and of cell lines HEK193, A549 and HepG2 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat. Dog