



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05199 - Goat Anti-EGR2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: EGR2, early growth response 2 (Krox-20 homolog, Drosophila), KROX20, Krox-20 homolog, KROX-20, Drosophila, homolog (early growth response-2), early growth response 2 (Krox-20 (Drosophila) homolog)

Official Symbol: EGR2

Accession Number(s): NP_000390

Human GeneID(s): [1959](#)

Immunogen

Peptide with sequence GGPLAPCSSRTRTP, from the C Terminus of the protein sequence according to NP_000390.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Weak band observed in 293 at approx 50kDa with 1 hour incubation of primary.

Species Reactivity

Tested:

Expected from sequence similarity: Human