



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06751-T - Goat Anti-MAPK3/ERK1 Antibody - Trial

Size: 20µg specific antibody in 40µl



Target Protein

Principal Names: MAPK3, ERK1, mitogen-activated protein kinase 3, HGNC:6877, P44ERK1, P44MAPK, PRKM3, protein kinase, mitogen-activated 3 (MAP kinase 3 p44), HS44KDAP, HUMKER1A, MGC20180, OTTHUMP00000174538, extracellular signal-regulated kinase 1, extracellular signal-related kinase 1

Official Symbol: MAPK3

Accession Number(s): NP_002737.2; NP_001035145.1; NP_001103361.1

Human GeneID(s): [5595](#)

Important Comments: No cross-reactivity expected with ERK2

Immunogen

Peptide with sequence GGEPRRTEGVGP-C, from the N Terminus of the protein sequence according to NP_002737.2; NP_001035145.1; NP_001103361.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Approx 40kDa band observed in Human Cerebellum lysates and in lysates of cell lines A431 and HepG2, and approx 42kDa observed in Human Frontal Cortex lysates (calculated MW of 40.1kDa according to NP_001035145.1 and 43.1kDa according to NP_002737.2). Recommended concentration: 0.1-1µg/ml. Primary incubation 1 hour at room temperature.


Immunofluorescence: Strong expression of the protein seen in the cytoplasm of HeLa and A431 cells. Recommended concentration: 10µg/ml.

Flow Cytometry: Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.


Species Reactivity

Tested: Human


Expected from sequence similarity: Human




EB06751 (0.1µg/ml) staining of Human Cerebellum (A) and Frontal Cortex (B) lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.




EB06751 (1µg/ml) staining of A431 (A) and (0.3µg/ml) HepG2 (B) cell lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.



EB06751 Immunofluorescence analysis of paraformaldehyde fixed HeLa cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing cytoplasmic and nuclear staining. Actin filaments were stained with phalloidin (red) and the nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB06751 Immunofluorescence analysis of paraformaldehyde fixed A431 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing cytoplasmic staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).



EB06751 Flow cytometric analysis of paraformaldehyde fixed A431 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control: Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.