



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08559 - Goat Anti-ERN1 / IRE1a Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ERN1, IRE1, endoplasmic reticulum to nucleus signaling 1, FLJ30999, IRE1P, MGC163277, MGC163279, ER to nucleus signalling 1, endoplasmic reticulum to nucleus signalling 1, inositol-requiring 1, protein kinase/endoribonuclease, IRE1a, IRE1-alpha

Official Symbol: ERN1

Accession Number(s): NP_001424.3

Human GeneID(s): [2081](#)

Non-Human GeneID(s): 78943 (mouse), 498013 (rat)

Immunogen

Peptide with sequence C-PNAHGKIKAMISD, from the internal region of the protein sequence according to NP_001424.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 110kDa band observed in Human Pancreas and Duodenum lysates (calculated MW of 110kDa according to NP_001424.3). Recommended concentration: 0.1-0.3µg/ml. Primary incubation 1 hour at room temperature.

Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB08559 (0.3 μ g/ml) staining of Human Duodenum (A) and (0.1 μ g/ml) Pancreas (B) lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence



EB08559 Flow cytometric analysis of paraformaldehyde fixed HeLa cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10 μ g/ml) followed by Alexa Fluor 488 secondary antibody (1 μ g/ml). IgG control: Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.