



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07568 - Goat Anti-FADS1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FADS1, fatty acid desaturase 1, D5D, FADS6, FADSD5, FLJ38956, FLJ90273, LLCDL1, TU12, delta-5 desaturase, delta-5 fatty acid desaturase, linoleoyl-CoA desaturase (delta-6-desaturase)-like 1

Official Symbol: FADS1

Accession Number(s): NP_037534.3

Human GeneID(s): [3992](#)

Non-Human GeneID(s): 76267 (mouse), 84575 (rat)

Immunogen

Peptide with sequence C-QPSFEPTKNKELTDE, from the internal region of the protein sequence according to NP_037534.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 50kDa band observed in Human Lung lysates (calculated MW of 51.9kDa according to NP_037534.2). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB07568 (0.3 μ g/ml) staining of Human Lung lysate (35 μ g protein in RIPA buffer). Primary incubation was 1 hour.
Detected by chemiluminescence.