

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06048 - Goat Anti-FBXL10 / PCCX2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: KDM2B, lysine (K)-specific demethylase 2B, FBL10, FBXL10, F-box and leucine-rich repeat protein 10, CXXC2, Fbl10, JHDM1B, PCCX2, JEMMA (Jumonji domain, EMSY-interactor, methyltransferase motif) protein, protein containing CXXC

domain 2, jumonji C domain-containing histone demethylase 1B

Official Symbol: KDM2B

Accession Number(s): NP_115979.3

Human GeneID(s): 84678

Non-Human GenelD(s): 30841 (mouse), 304495 (rat)

Important Comments: This antibody is expected to recognize isoform a (NP_115979.3)

only.

Immunogen

Peptide with sequence C-FGQVEEKLLQKLS, from the C Terminus of the protein sequence according to NP_115979.3.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:32000.

IHC: Paraffin embedded Human Skin. Recommended concentration: 5-7µg/ml.


Species Reactivity

Tested: Human


Expected from sequence similarity: Human, Mouse, Rat, Dog



EB06048 ($5\mu g/ml$) staining of paraffin embedded Human Skin. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



EB06048 ($7\mu g/ml$) staining of paraffin embedded Human Skin. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB06048 Negative Control showing staining of paraffin embedded Human Skin, with no primary antibody.