



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06331 - Goat Anti-FEZ2 / Zygin 2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FEZ2, HUM3CL, fasciculation and elongation protein zeta 2 (zygin II), pre-T/NK cell associated protein (3CI), fasciculation and elongation protein zeta 2, MGC117372, zygin 2

Official Symbol: FEZ2

Accession Number(s): NP_005093.1

Human GeneID(s): [9637](#)

Immunogen

Peptide with sequence GERRGRDGLGAPT-C, from the N Terminus of the protein sequence according to NP_005093.1.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Western Blot: Preliminary experiments gave no signal but low background in Human Brain and Lung extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human