

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB07743 - Goat Anti-Flotillin 2 / FLOT2 (internal) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FLOT2, flotillin 2, reggie-1, ESA, ESA1, M17S1, Flotillin 2 (epidermal

surface antigen 1)

Official Symbol: FLOT2

Accession Number(s): NP_004466.2

Human GenelD(s): 2319

Non-Human GenelD(s): 14252 (mouse), 83764 (rat)

Immunogen

Peptide with sequence C-QIYQDRDQFAK, from the internal region of the protein sequence according to NP_004466.2.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Western Blot: Approx 37kDa band observed in Human, Mouse and Rat Brain lysates (calculated MW of 41.7kDa according to NP_004466.1). Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog