



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08625 - Goat Anti-FPRL1 (aa13-26) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FPRL1, formyl peptide receptor-like 1, ALXR, FMLP-R-II, FMLPX, FPR2A, FPRH1, FPRH2, HM63, LXA4R, lipoxin A4 receptor (formyl peptide receptor related)

Official Symbol: FPRL1

Accession Number(s): NP_001453.1

Human GeneID(s): [2358](#)

Important Comments: Reported variants represent identical protein: NP_001005738.1 and NP_001453.1.

Immunogen

Peptide with sequence EEVSYESAGYTVLR-C, from the N Terminus of the protein sequence according to NP_001453.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary testing showed a band at approx 40kDa in Human Heart and approx. 37-38kDa in Human Spleen lysate after 1µg/ml antibody staining (calculated MW of 39.0kDa according to NP_001453.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human