

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08625 - Goat Anti-FPRL1 (aa13-26) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FPRL1, formyl peptide receptor-like 1, ALXR, FMLP-R-II, FMLPX, FPR2A, FPRH1, FPRH2, HM63, LXA4R, lipoxin A4 receptor (formyl peptide receptor

related

Official Symbol: FPRL1

Accession Number(s): NP_001453.1

Human GeneID(s): 2358

Important Comments: Reported variants represent identical protein: NP_001005738.1

and NP_001453.1.

Immunogen

Peptide with sequence EEVSYESAGYTVLR-C, from the N Terminus of the protein sequence according to NP_001453.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Preliminary testing showed a band at approx 40kDa in Human Heart and approx. 37-38kDa in Human Spleen lysate after 1µg/ml antibody staining (calculated MW of 39.0kDa according to NP_001453.1). Primary incubation 1 hour at room temperature.

Species Reactivity

Tested:

Expected from sequence similarity: Human