

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB05205 - Goat Anti-FYB Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: FYB, FYN binding protein (FYB-120/130), ADAP, SLAP-130,

FYN-binding protein (FYB-120/130)

Official Symbol: FYB

Accession Number(s): NP_001456

Human GeneID(s): 2533

Immunogen

Peptide with sequence YDDIADGCIYDND, from the C Terminus of the protein sequence according to NP_001456.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: No signal obtained yet but low background observed in Jurkat and U937 lysates at upto 2μg/ml.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse