



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09856 - Goat Anti-Glutamate Dehydrogenase Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GDH, GDH1, GLUD, glutamate dehydrogenase 1, MGC132003, GLUD1

Official Symbol: GLUD1 / GLUD2

Accession Number(s): NP_005262.1; NP_036216.2

Human GeneID(s): [2746](#), [2747](#)

Non-Human GeneID(s): 14661 (mouse), 24399 (rat)

Important Comments: This antibody is expected to recognize both reported isoforms (NP_005262.1; NP_036216.2).

Immunogen

Peptide with sequence C-ESEEQKRNVRGILR, from the internal region of the protein sequence according to NP_005262.1; NP_036216.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 55kDa band observed in Human Brain (Cerebellum) and Mouse Brain lysates and in Human, Mouse and Rat Liver lysates (calculated MW of 61.4kDa according to Human NP_005262.1, of 61.3kDa according to Mouse NP_032159.1 and of 61.4kDa according to Rat NP_036702.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.

IHC: In paraffin embedded Human Cerebral Cortex shows mitochondrial staining in neuronal cell bodies and axons. Recommended concentration: 3-5µg/ml.

Immunofluorescence: Customer found particulate cytoplasm staining in HeLa.

Additional validation: This antibody has been successfully used in the following paper: Sikorski et al. (2018) PMID: 30377371.

Species Reactivity


Tested: Human, Mouse, Rat

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow


Specific Reference

This antibody has been successfully used in the following paper:


Krzysztof Sikorski, Adi Mehta, Marit Inngjerdigen, Flourina Thakor, Simon Kling, Tomas Kalina, Tuula A. Nyman, Maria Ekman Stensland, Wei Zhou, Gustavo A. De Souza, Lars Holden, Jan Stuchly, Markus Templin and Fridtjof Lund-Johansen
A high-throughput pipeline for validation of antibodies
Nat Methods. 2018 Nov;15(11):909-912
PMID: 30377371



EB09856 (0.03 μ g/ml) staining of Human (A) and Mouse (B) Brain lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB09856 (0.01 μ g/ml) staining of Human (A), Mouse (B) and Rat (C) Liver lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB09856 (3.75µg/ml) staining of paraffin embedded Human Cerebral Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.