



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB07187 - Goat Anti-GNB3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GNB3, guanine nucleotide binding protein (G protein), beta polypeptide 3, G protein, beta-3 subunit, GTP-binding regulatory protein beta-3 chain, guanine nucleotide-binding protein G(I)/G(S)/G(T) beta subunit 3, guanine nucleotide-binding protein, beta-3 subunit, hypertension associated protein, transducin beta chain 3

Official Symbol: GNB3

Accession Number(s): NP_002066.1

Human GeneID(s): [2784](#)

Important Comments: No cross-reactivity expected with GNB1, 2 and 4

Immunogen

Peptide with sequence DCMSLAVSPDFN, from the internal region of the protein sequence according to NP_002066.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Approx 36kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 37.2kDa according to NP_002066.1). Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB07187 (1.0µg/ml) staining of Human Brain (Frontal Cortex) lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.