

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06142 - Goat Anti-GRAF / ARHGAP26 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ARHGAP26, GRAF, Rho GTPase activating protein 26, KIAA0621, FLJ42530, OPHN1L, OPHN1L1, GTPase regulator associated with FAK, GTPase regulator associated with focal adhesion kinase pp125(FAK), GTPase regulator

associated with the focal adhesion kinase pp125 **Official Symbol:** ARHGAP26

Accession Number(s): NP_055886.1; NP_001129080.1

Human GeneID(s): 23092

Non-Human GenelD(s): 71302 (mouse)

Important Comments: This antibody is expected to recognise both reported isoforms.

Immunogen

Peptide with sequence C-KTGLIPENYVEFL, from the C Terminus of the protein sequence according to NP_055886.1; NP_001129080.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:64000.

Western blot: Approx 90kDa band observed in Human Heart lysates (calculated MW of 92.2kDa according to NP_055886.1). Recommended concentration: 0.3-1μg/ml. An additional strong band of unknown identity was also consistently observed at 16kDa. This band was successfully blocked by incubation with the immunising peptide.

IHC: Paraffin embedded Human Colon. Recommended concentration: 3.75µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB06142 ($0.3\mu g/ml$) staining of Human Heart lysate ($35\mu g$ protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Primary incubation was 1 hour. Detected by chemiluminescence.



EB06142 (3.75μg/ml) staining of paraffin embedded Human Colon. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.