



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05080 - Goat Anti-GRB14 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: GRB14, growth factor receptor-bound protein 14

Official Symbol: GRB14

Accession Number(s): NP_004481

Human GeneID(s): [2888](#)

Immunogen

Peptide with sequence KLKHYCARIAL, from the C Terminus of the protein sequence according to NP_004481.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western blot: Approx 50-55kDa, 40kDa and 30kDa bands seen in HeLa cell lysates. Recommended for use at 0.1-0.5µg/ml.

has to fail as unexplainable 40 and 30kDa bands

Species Reactivity

Tested: Human

Expected from sequence similarity: Human